# PHILIPS

# RESPIRONICS

# Amara View minimal contact full-face mask Instructions for Use

# Intended Use

This mask is intended to provide an interface for application of CPAP or bi-level therapy to patients. The mask is for single patient use in the home or multi-patient use in the hospital/institutional environment. The mask is to be used on patients (>66lbs/30kg) for whom CPAP and bi-level therapy has been prescribed.

Caution: U.S. Federal law restricts this device to sale by or on the order of a physician.

# Note<sup>-</sup>

- An exhalation port is built into this mask so a separate exhalation port is not required.
- This mask is not made with natural rubber latex or DEHP.

# Symbols Glossary

Refer to http://www.symbols.philips.com for a description of the symbols used on this device and its packaging.









Not Made with Natural

🚱 Tip

X1 System One Resistance Control

# A Warnings


- This mask is not suitable for providing life support ventilation.
- This mask is designed for use with CPAP or bi-level systems recommended by your health care professional or respiratory therapist. Do not wear this mask unless the CPAP or bi-level system is turned on and operating properly. Do not block or try to seal the exhalation port. Explanation of the Warning: CPAP systems are intended to be used with special masks with connectors which have vent holes to allow continuous flow of air out of the mask. When the CPAP machine is turned on and functioning properly, new air from the CPAP machine flushes the exhaled air out through the attached mask exhalation port. However, when the CPAP machine is not operating, enough fresh air will not be provided through the mask, and exhaled air may be rebreathed.


Rebreathing of exhaled air for longer than several minutes can in some circumstances lead to suffocation. This warning applies to most models of CPAP systems.

- If oxygen is used with the device, the oxygen flow must be turned off when the device is not operating. Explanation of the Warning: When the device is not in operation, and the oxygen flow is left on, oxygen delivered into the ventilator tubing may accumulate within the device enclosure. Oxygen accumulated in the device enclosure will create a risk of fire.
- At a fixed flow rate of supplemental oxygen flow, the inhaled oxygen concentration will vary, depending on the pressure settings, patient breathing pattern, mask selection, and the leak rate. This warning applies to most types of CPAP and bi-level machines.
- At low CPAP or EPAP pressures the flow through the exhalation port may be inadequate to clear all exhaled gas from the tubing. Some rebreathing may occur.
- Some users may experience skin redness, irritation, or discomfort. If this happens, discontinue use and contact your healthcare professional.
- Consult a physician if any of the following symptoms occur: Unusual chest discomfort, shortness of breath, stomach distension, belching, severe headache, blurred vision, drying of the eyes, eye pain or eye infections.
- This mask should not be used on patients who are uncooperative, obtunded, unresponsive, or unable to remove the mask.
- This mask is not recommended if the patient is taking a prescription drug that may cause vomiting.
- Consult a physician or dentist if you encounter tooth, gum, or jaw soreness. Use of a mask may aggravate an existing dental condition.
- A minimum of  $3 \text{cm H}_2\text{O}$  (hPa) must be maintained when using this mask.
- Attaching an exhalation device requires therapy pressure level adjustment to compensate for increased leak.
- Do not overtighten the headgear straps. Watch for signs of overtightening, such as
  excessive redness, sores or bulging skin around the edges of the mask. Loosen the
  headgear straps to alleviate symptoms.
- Do not block or seal the anti-asphyxia valve or exhalation ports.

# If your mask uses magnetic headgear clips, the following warnings apply:

- The mask assembly contains magnets. Some medical devices may be affected by magnetic fields. The magnetic clips in this mask should be kept at least 2 in. (50 mm) away from any active medical device with special attention to implanted devices such as pacemakers, defibrillators and cochlear implants.
- Do not use in or near magnetic resonance imaging (MRI) equipment.
- Keep unassembled magnetic headgear clips out of reach of children. In case of accidental swallowing, seek medical assistance immediately.





# Figure (1)

- (A) Mask Cushion with Exhalation Ports (Do not block)
- (B) Mask Frame
- č Elbow with Anti-Asphyxia Valve (Do not block)
- Ouick Release Tube with swivel
- Headgear Crown Strap
- D Quick Release Tabs
   Quick Release Tube
   F Headgear Crown Sti
   G Headgear Top Side 2
   H Headgear Bottom S
   Mask Frame Top He Headgear Top Side Strap
- Headgear Bottom Side Strap
- Mask Frame Top Headgear Slot
- Magnetic Headgear Clip Receiver
- Magnetic Headgear Clip


# **Optional:**

L Talon Headgear Clip

# Verify the Anti-Asphyxia Valve

The anti-asphyxia valve consists of an air inlet and a flapper (2). With the airflow turned off, verify that the elbow with the flapper is lying flat (2) so that room air can flow in and out through the air inlet. Next, with the airflow on, the flapper should now cover the air inlet and air from the CPAP or bi-level device should flow into the mask (2). If the flapper does not close or does not function properly, replace the mask.

▲Warning: Do not block or seal the anti-asphyxia valve.



# Contraindications

This mask may not be suitable for use on patients with the following conditions: recent eye surgery or dry eyes, hiatal hernia, excessive reflux, impaired cough reflex, impaired cardiac sphincter function; or on patients unable to remove the mask by themselves.

### Before Use Read and Understand the Instructions Completely.

- Hand wash the mask.
- Wash your face.

**A** Caution: Do not use moisturizer/lotion on your hands or face.

- Inspect the mask and replace if the cushion has hardened or is torn, or if any parts are broken.
- Verify that the therapy device, i.e., ventilator, including the alarms and safety systems, has been validated prior to use.
- Verify therapy device pressure(s).

# **Cleaning Instructions**

Hand wash the mask before first use and daily. The headgear should be hand washed weekly. The headgear does not need to be removed for daily cleaning of the mask.

1. Hand wash mask and headgear in warm water with liquid dishwashing detergent.

2. Rinse thoroughly. Air dry completely before use. Make sure the mask is dry before use. Lay the headgear flat or line dry.

Caution: Do not place the headgear into the dryer.

**Caution:** Do not use bleach, alcohol, cleaning solutions containing bleach or alcohol, or cleaning solutions containing conditioners or moisturizers.

Caution: Any deviation from these instructions may damage the product.


**Caution:** Inspect the mask for damage or wear (cracking, crazing, tears, etc). Discard and replace any components as necessary.

#### Institutional Disinfection

Clean and disinfect the mask daily and between patients when in use in a clinical setting. Refer to the Disinfection Guide for Professional Users to reprocess between patients in a clinical setting. Access the latest version of the Disinfection Guide at www.healthcare. philips.com or by contacting customer service at 1-800-345-6443 (USA or Canada only) or +1-724-387-4000.

# Achieving the Right Fit Sizing the Mask

The mask should rest comfortably against your face. The cushion nostril opening goes under your nostrils (3). The headgear should be away from your eyes. The bottom of the mask cushion should rest just above your chin with your mouth slightly open.




# Before putting on the mask:

1. Headgear - Loosen to a large setting.

2. Magnetic Headgear Clips - Grip and twist outward to disconnect from the magnetic headgear clip receiver. ④

 ${\rm Optional:}$  If your mask uses talon headgear clips, grip and twist away from the frame to disconnect.  $\textcircled{\begin{tabular}{ll} \end{tabular}}$ 



# Adjusting the mask:

1. Hold the mask against your face. Pull the headgear over your head (5).


2. Mask cushion - Place your nose over the cushion nostril opening. Push the mask cushion up.

Solution of the second second

3. Magnetic Headgear Clips - Connect the magnetic headgear clip to the magnetic headgear clip receiver.

**Optional:** If your mask uses talon headgear clips, press into place to connect. 4. Headgear - Pull back the tabs to adjust the top and bottom side straps. The headgear should fit loosely against your face. Do not overtighten. **6**7

5. Headgear - Pull the crown strap tab over to reduce leaks around your nose (8).



6. Headgear - Make final adjustments by pulling back the tabs to evenly adjust the top and bottom side straps. The headgear should feel snug, but not uncomfortably tight (6)(7).

**Note:** You should be able to fit one of your fingers underneath the top and bottom side straps.

# Using the mask:

1. Connect the flexible tubing (included with your CPAP or bi-level device) to the quick release tube.

2. Lie down. Turn the therapy device on. Breathe normally. Keep your mouth slightly open.

3. Assume different sleeping positions. Move around until comfortable. If there are any excessive air leaks, make final adjustments. If there are any excessive air leaks around the nose, adjust the crown strap first. This will move the mask closer or further away from the bottom of your nose. Some air leaking is normal.

# 🔂 Comfort Tips

- The most common mistake is overtightening the headgear. The headgear should fit loose and comfortable. If your skin bulges around the mask or if you see red marks on your face, loosen the headgear.
- Adjust the crown headgear strap to reduce leaks around your nose and to lift the straps off of your ears (8).
- Adjust the top and bottom side straps to reduce leaks at the sides or bottom of the mask cushion (6)(7).

### Removing the Mask

Disconnect the magnetic headgear clips and pull over the top of your head.

**Optional:** If your mask uses talon headgear clips, disconnect the headgear clips and pull over the top of your head.

### Disassembly


1. Magnetic Headgear Clips - Grip and twist away from the magnetic headgear clip receiver to disconnect.

**Optional:** If your mask uses talon headgear clips, grip and twist away from the frame to disconnect.


2. Headgear - Undo the straps. Pull out of headgear clips and mask frame top headgear slots.

3. Mask Cushion - Grasp the top and the bottom of the mask cushion, and pull it off the mask frame 9.

**A** Caution: Do not pull the mask cushion by the nostril cushion opening.



4. Quick Release Tube - Squeeze the quick release tabs. Pull the tube away from the mask frame 10.



#### Assembly

1. Headgear - Thread the headgear straps through the mask frame top headgear slots and magnetic headgear clips.

2. Magnetic Headgear Clips - Connect the headgear clip to the magnetic headgear clip receiver.

**Optional:** If your mask uses talon headgear clips, push the headgear clips onto the mask frame.

3. Mask Cushion - Grasp both sides and push it onto the mask frame.

4. Quick Release Tube - Push onto the elbow until the quick release tabs click into place.


### Philips Respironics System One Resistance Control

Your mask when combined with a Philips Respironics System One device, provides optimal resistance compensation. This mask value is X1 and should be set by your provider. Notes:

- Compare the mask to the device. See your device manual if the values do not match.
- System One is not compatible with masks that use a separate exhalation device.

# Specifications

▲ Warning: The technical specifications of the mask are provided for your healthcare professional to determine if it is compatible with your CPAP or bi-level therapy device. If used outside these specifications, or if used with incompatible devices, the mask may be uncomfortable, the seal of the mask may not be effective, optimum therapy may not be achieved, disconnect alarms may not be detected, and leak, or variation in the rate of leak, may affect device function. The pressure flow curve shown below is an approximation of expected performance. Exact measurements may vary.



#### Pressure Flow Curve

Resistance with Anti-Asphyxia Valve Closed to Atmosphere Drop in Pressure at

	50 SLPM	100 SLPM
All sizes	$0.8 \mathrm{cm}\mathrm{H_2O}$	$2.0 \text{ cm H}_2\text{O}$

**Closed to Atmosphere:** A minimum of 3 cm  $H_2O$  must be maintained to ensure that the anti-asphyxia valve closes.

**Open to Atmosphere:** In the absence of positive airway pressure, the antiasphyxia valve will open.

# Inspiratory and Expiratory Resistance with Anti-Asphyxia Valve Open to Atmosphere at 50 SLPM

Inspiratory Resistance	
Expiratory Resistance	

1.0 cm H<sub>2</sub>O 1.5 cm H<sub>2</sub>O

#### Deadspace

s .	150.4 mL
М	148.8 mL
L	167.8 mL

# Sound Levels

A-weighted Sound Power Level	32.5 dBA
A-weighted Sound Pressure Level at 1m	24.2 dBA

## Disposal

Dispose of in accordance with local regulations.

Dispose of mask and headgear annually.

### Storage Conditions

Temperature: -4° F to 140° F (-20° C to 60° C) Relative Humidity: 15% to 95%, non-condensing

### LIMITED WARRANTY

Respironics, Inc., a Philips company, warrants that the system shall be free from defects of workmanship and materials and will perform in accordance with the product specifications for a period of ninety (90) days from the date of sale by Respironics, Inc. to the dealer. If the product fails to perform in accordance with the product specifications, Respironics, Inc. will pay customary freight charges from Respironics, Inc. to the dealer location only. This warranty does not cover damage caused by accident, misuse, abuse, alteration, water ingress, and other defects not related to material or workmanship. The Respironics, Inc. Service department shall examine any devices returned for service, and Respironics, Inc. reserves the right to charge an evaluation fee for any returned device as to which no problem is found after investigation by Respironics, Inc. Service.

This warranty is non-transferable by unauthorized distributors of Respironics, Inc. products and Respironics, Inc. reserves the right to charge dealers for warranty service of failed product not purchased directly from Respironics, Inc. or authorized distributors.

Respironics, Inc. disclaims all liability for economic loss, loss of profits, overhead, or consequential damages which may be claimed to arise from any sale or use of this product. Some states do not allow the exclusion or limitation of incidental or consequential damages, so the above limitation or exclusion may not apply to you.

This warranty is given in lieu of all other express warranties. In addition, any implied warranties – including any warranty of merchantability or fitness for the particular purpose – are limited to two years. Some states do not allow limitations on how long an implied warranty lasts, so the above limitation may not apply to you. This warranty gives you specific legal rights, and you may also have other rights which vary from state to state. To exercise your rights under this warranty, contact your local authorized Respironics, Inc. dealer or contact Respironics, Inc. at:

1001 Murry Ridge Lane

Murrysville, Pennsylvania 15668-8550

1-724-387-4000



Respironics Inc. 1001 Murry Ridge Lane Murrysville, PA 15668 USA

