

A man and a woman are sleeping peacefully in a bed. The man, in the foreground, is wearing glasses and a clear nasal CPAP mask connected to a white tube. The woman is sleeping behind him, her face partially visible. The background is a soft, out-of-focus light blue wall.

PHILIPS

RESPIRONICS

Wisp

Minimal contact nasal mask

Performance has a
new style

Wisp innovation

Efficient and simple

Developing solutions that meet the needs of a rapidly changing industry calls for fresh thinking and collaboration. At Philips Respiroics, we have traveled around the world asking clinicians, homecare providers, and patients how we could improve their experience.

Clinicians wanted us to resolve the tradeoff between nasal and pillows masks. Homecare providers wanted a solution that could reduce operational complexity by getting more patients served from one package. And patients wanted a better overall mask experience that improved their lifestyle.

When you put it all together, the answer is Wisp.

Traditional nasal mask

- Larger footprint
- Intimidation factor

- + Seal/stability
- + Comfort

Wisp

- + Minimal size
- + Lightweight

- Pressure restrictions
- Nasal discomfort

Traditional pillows mask





Magnetic personality

Magnetic clip accessory for enhanced ease of use

The easy-to-use Wisp magnetic clip is standard in the Wisp fitpack and an accessory part that can be retrofitted to existing Wisp masks. It can replace the standard headgear clip to enhance ease of use for patients with limited dexterity.






Wisp features

Self-adjusting top strap
with tube holder loop

Reversible fabric frame of soft
suedette and silky sateen (also
available in a clear silicone)

Magnetic clips for
easy adjustment





Tip-of-the-nose
cushion design

Quick-release
elbow tabs

Wisp

Wisp patient fitting guide

Follow the instructions to get started with fitting the Wisp nasal mask.

- 

1 Before putting on the mask, select the cushion that fits the width of the nose without blocking the nostrils. The top of the cushion may sit between the bridge and tip of the nose.

Note: Use the mask sizing gauge to select the cushion that fits the width of the nose without blocking the nostrils.
- 

2 Disconnect the headgear tab or clip. Peel the headgear tabs away from the headgear strap and pull them through the slot in the mask frame, or if using the clip, twist outward.
- 

3 Pull the headgear over the head and place the cushion over the nose.

Note: The headgear should lie flat against the back of the head.
- 

4 Thread the tab through the slot in the mask frame or press the clip against the headgear slot to reattach.

Note: If the clips are not used, the headgear tabs can be threaded through the slots on the mask frame.
- 

5 Peel the headgear tabs away from the headgear straps and then press them back together to reattach. Connect the side straps to the headgear tab located at the base of the crown strap. Connect the bottom strap tabs near the back of the headgear. Adjust the side and bottom straps evenly.

Note: Do not over-tighten. The headgear should feel loose and comfortable.
- 

6 Connect the Wisp tubing and elbow swivel to the front of the mask. Optionally, slide the mask tubing through the headgear crown loop before pulling on the headgear.
- 

7 Assume different sleeping positions. Move around until comfortable. If there are any excessive air leaks, make final adjustments while lying down.

Note: Using the headgear crown loop can help make some sleeping positions much more comfortable.

Choose which Wisp works

Wisp features a unique, compact design that is created for maximum comfort. The “**tip-of-the-nose**” cushion is designed to sit below the nose bridge, making minimal contact with the face while creating an optimal seal. Wisp’s sleek design eliminates the need for a forehead pad, offering a better field of vision.



Petite



Small/medium



Large



Extra large

We’ve taken the guesswork out of finding the fit that’s right. The fitpack includes small/medium, large, and extra large sizes, with the petite cushion being sold separately.

Choice of frame




Clear frame

The silicone frame is lightweight and flexible, offering a smooth feel for patients who prefer the look of clear silicone.

Fabric frame

The reversible fabric frame of soft suedette and silky sateen allows patients to choose which side feels best.

DreamMapper



DreamMapper is the #1 downloaded PAP therapy engagement app¹

When it comes to your patients reaching their sleep therapy goals, DreamMapper² is there for them with the tools, knowledge and resources they can use to stay motivated.

Visit: www.philips.com/mydreammapper to experience DreamMapper.

¹Source GooglePlay, 5/9/2016.

²To see which therapy devices are compatible with DreamMapper, visit: www.dreammapper.com/compatible.

Ordering information

Wisp fitpack with headgear; small/medium, large, extra large cushions

| | |
|---|---------|
| Wisp mask, clear frame, headgear, magnetic clips | 1094050 |
| Wisp mask, fabric frame, headgear, magnetic clips | 1094051 |

Wisp fitpack without headgear; small/medium, large, extra large cushions

| | |
|---|---------|
| Wisp mask, clear frame, without headgear | 1094055 |
| Wisp mask, fabric frame, without headgear | 1094056 |

Wisp mask, clear frame without headgear

| Part number | Part number |
|--------------|-------------|
| Petite | 1101550 |
| Small/medium | 1101551 |
| Large | 1101552 |
| Extra large | 1113957 |

Wisp mask, fabric frame without headgear

| Part number | Part number |
|--------------|-------------|
| Petite | 1101553 |
| Small/medium | 1101554 |
| Large | 1101555 |
| Extra large | 1113958 |

Wisp mask, clear frame with headgear

| Part number | Part number |
|--------------|-------------|
| Petite | 1118060 |
| Small/medium | 1118061 |
| Large | 1118062 |
| Extra large | 1118063 |

Wisp mask, fabric frame with headgear

| Part number | Part number |
|--------------|-------------|
| Petite | 1118064 |
| Small/medium | 1118065 |
| Large | 1118066 |
| Extra large | 1118067 |

Wisp youth fitpack with reduced-size headgear; small/medium, large, extra large cushions

| | |
|---|---------|
| Wisp mask, fabric frame, headgear, magnetic clips | 1109298 |
|---|---------|

Replacement parts

| Part number | Part number |
|---------------------------|-------------|
| Petite Wisp cushion | 1094086 |
| Small/medium Wisp cushion | 1094087 |
| Large Wisp cushion | 1094088 |
| Extra large Wisp cushion | 1112031 |



Replacement parts

| Part number | Part number |
|-----------------------------|-------------|
| Reduced-size Wisp headgear* | 1109307 |
| Large Wisp headgear* | 1094078 |
| Wisp headgear* | 1094082 |



*Clips and cushions sold separately

Replacement parts

| Part number | Part number |
|------------------------------|-------------|
| Wisp headgear clips | 1094083 |
| Wisp magnetic headgear clips | 1120434 |



Replacement parts

| Part number | Part number |
|--------------------|-------------|
| Wisp clear frame* | 1094080 |
| Wisp fabric frame* | 1094081 |



*Cushion sold separately

Learn more

For more information about Wisp, visit www.philips.com/wisp.





© 2017 Koninklijke Philips N.V. All rights reserved.
Specifications are subject to change without notice.

www.philips.com/respironics

Caution: US federal law restricts these devices to
sale by or on the order of a physician.

RRDPGH EB 12/1/17 MCI 4108115 PN 1103630
1010 Murry Ridge Lane, Murrysville, PA 15668
800-345-6443 • 724-387-4000