



Invacare HomeFill Oxygen System

Pulse dose and continuous flow options available

Technical Specifications

ML9 Tank

Size: Lasts up to 10 hours on 2lpm/20bpm
Weight: 4lbs

"D" Tank

Size: Lasts 3 hours on 2lpm continuous
Features: Integrated Regulator

- Lightweight**
- Easy to Use**
- Easy to Fill**



Respironics SimplyGo

The Respironics SimplyGo portable oxygen concentrator features continuous flow and pulse-dose delivery in a sleek 10-pound package. For the size, the SimplyGo has more than twice the oxygen output of other small POC options on the market.

- Prescription is required for this product
- Covered by most insurances

Technical Specifications

Size:	29.2 cm x 25.4 cm x 15.2 cm (with standard battery)
Weight:	10 lbs. (4.5 kg) with battery installed
User Interface:	Touch pad
Operating temperature:	-4° F to 140° F unit only
Duration:	3 hours (pulse dose setting of 2 LPM) 0.9 hours (continuous flow setting of 2 LPM)
Demand Continuous flow:	Up to a setting of 2 LPM
Demand Pulse flow:	Up to a setting of 6 LPM



Respironics SimplyGo Mini

The Respironics SimplyGo Mini portable oxygen concentrator is the smallest and lightest portable oxygen concentrator offered at Home Oxygen Company. Called the most advanced portable oxygen concentrator ever made, the SimplyGo Mini offers color LED touch display, easy access to flow settings and weighs only 5 pounds.

- Prescription is required for this product
- Covered by most insurances

Technical Specifications

Size:	23.9 cm x 21.1 cm x 9.1 cm (with standard battery)
Weight:	5 lbs. (2.3 kg) with standard battery installed 6 lbs. (2.7 kg) with extended battery installed
User Interface:	Touch screen LCD
Operating temperature:	41° F to 95° F
Duration:	Standard battery lasts up to 4.5 hours (pulse setting of 2, 20 BPM) Extended battery lasts up to 9 hours (pulse setting of 2, 20 BPM)