

drive

 DeVilbiss[®]
HEALTHCARE



iGo[®]2

DISCOVER NEW HORIZONS

DISCOVER NEW HORIZONS IN PORTABLE OXYGEN THERAPY

Through its groundbreaking SmartDose technology, the iGo2 has created a completely new class within portable oxygen concentrators. The iGo2 automatically adjusts the oxygen delivery bolus rate to the patient's changing breathing frequency.

Exclusive SmartDose Technology Optimising Oxygen Delivery

The **patented SmartDose technology** enables breath-rate-controlled oxygen therapy. When in SmartDose mode, the algorithm monitors patient breath rate changes, therefore manual adjustment is no longer necessary, thus promoting an active and individual lifestyle. Oxygen is automatically adjusted and increased based on user activity. During periods of rest, the device will revert to the baseline setting.



SMARTDOSE PROMOTES
AN ACTIVE AND
INDIVIDUAL LIFESTYLE.



EFFICIENT THERAPY - OPTIMUM DEMAND



SmartDose™ Technology

Automatically adjusts oxygen production & delivery in real time to meet patient demand.



Precise Breath Detection

Due to its advanced trigger technology, the iGo2 accurately detects the patient's breathing.



Lightweight & Compact

One of the smallest and lightest portable oxygen concentrators on the market.



Battery Exchange

The iGo2 battery pack can be swapped within seconds, even whilst the unit is inside the carry case.

iGo[®]2

ALL AT A GLANCE.



The high-resolution colour display shows all essential information.

The user interface has been designed with an intuitive menu and easy to use control buttons to ensure user-friendly operation.

Adjusting the iGo2 settings to your preferences (+/- buttons) is effortless, as all controls are clearly identified.

Specified selection buttons for the SmartDose mode (SD) and alarm muting ensures safe and easy usage.



Optional Bluetooth Function (Bluetooth models only)



Displays SmartDose mode activation



O2 concentration purity indicator



Settings (1-5)



Battery status indicator



Air Travel
Approved



User-friendly LCD



Durable,
High Quality

The iGo2 portable oxygen concentrator impresses with its advanced and innovative design. Whilst small and compact, the device is both highly durable and robust. The iGo2 housing is constructed from impact-resistant material for added protection.

iGo[®]2 BLUETOOTH








 UNLOCK NEW HORIZONS


ACCESIBLE THERAPY DATA.

The iGo2 model with optional Bluetooth allows access to all the essential device information on your smartphone, using the SmartLink O2 application.

iGo[®]2 APP SmartLinkO2

YOUR BONUS WITH THE IGO2 BLUETOOTH

-  Displays remaining battery life (hours and minutes)
-  Displays oxygen quality
-  Displays therapy settings
-  Displays operating and usage hours
-  Displays system information; firmware version, serial number, and device type
-  Optional remote function via the MySmartLink Cloud Server
-  Easy device remote diagnosis (Service Provider)



SmartLinkO2



UNLIMITED INDEPENDENCE

The iGo2 portable oxygen concentrator provides simple and independent oxygen therapy wherever you are. Complete with interchangeable battery, car auxiliary cable, and universal power adapter, the iGo2 can be used at anytime, anywhere.

The carry case shoulder bag also included comfortably transports the iGo2 even during use yet remains virtually invisible.

Travelling becomes easier. Whether by car, bus, train, ship, or aeroplane, the iGo2 will accompany you on all your journeys and expeditions. **With the iGo2 you can experience new horizons and unlimited travel.**



COMING SOON*

Back-friendly travelling with the iGo2 backpack

The iGo2 optional backpack accessory provides prime protection for the device. Like the carry case, the backpack is light weight and also designed to accommodate iGo2 operation.

Thanks to its soft back padding, the backpack can be comfortably worn for long periods. Three various sized accessory pockets allow for storing the universal power adapter, car adapter, spare battery, and extra personal items.

*Expected to be available from autumn / winter 2024



TECHNICAL SPECIFICATIONS

**3,5
HOURS***
BATTERY
LIFE

*Setting 2 (20BPM)

**2,2 KG
WEIGHT**

**94%
OXYGEN
CONCENTRATION**



INCLUDED IN DELIVERY



- 1 CARRY CASE
- 2 BATTERY
- 3 POWER PLUG*
- 4 POWER SUPPLY UNIT
- 5 CAR AUXILIARY CABLE

*country specific

YOUR DRIVE DEVILBISS DEALER:



Further information can be found at: www.igo2-poc.de



Drive Medical GmbH & Co. KG • Leutkircher Straße 44 • 88316 Isny/Allgäu
DeVilbiss Healthcare GmbH • Kamenzer Straße 3 • 68309 Mannheim
T +49 (0) 7562 9724 0 • info@drivemedical.de • www.drivedevilbiss.de

DeVilbiss® iGo®2 is a registered trademark of DeVilbiss Healthcare.
SmartDose™ is a registered trademark of Inovo.